






Company	3Shape	Medit	Carestream™ Dental
<b>Product</b>	<b>3Shape TRIOS®</b> 	<b>Medit i500 &amp; i700</b> 	<b>CS 3700</b> 
<b>Scanner File Output</b>	DCM and STL	Open (meditmesh, OBJ, PLY, STL.)	STL and PLY
<b>Scanner Type</b>	LED Pattern Projection	3D-in-motion video technology	Parallel Confocal
<b>Scanner Accuracy (in microns)</b>	6.9	Full-Arch: 21.0µm ± 1.48	<9
<b>Infection Control</b>	Autoclavable/cold sterile.	Autoclavable	Autoclavable/cold sterile.
<b>Tooth Surface Spray Required?</b>	No	No	No
<b>Per Scan Click Fee?</b>	No	No	No
<b>Subscription (data package and/or software updates)</b>	\$1800 to \$3200/yr starting in yr 2	No subscription fees	No subscription fees (software updates included in initial purchase)
<b>Lapsed Time (min) until digital model file is available</b>	0 to 10; Varies depending on office's internet connection	0 to 5; Varies Depending on office's internet connection	0 to 10; Varies depending on office's internet connection
<b>Hardware Warranty?</b>	\$750/yr starting year 2 (optional)	1 year included with purchase	\$2,000/yr starting year 2 (optional)
<b>In-House Training</b>	1 to 1.5 days	Training available through Medit's partner network. Online assistance available.	1 day
<b>Suggested Retail Price</b>	\$22,400-\$50,600	\$12,000	\$30,000 - \$36,000
<b>Distributors</b>	Henry Schein, Biolase, Nashville Dental, Patterson Dental, and Great Lakes Ortho	Direct, Henry Schein, Pearson, and CAD-Ray	Direct & Benco
<b>Ortho applications currently?</b>	SureSmile and Clear Correct	Invisalign, Clear Correct, SureSmile	SureSmile and Clear Correct
<b>Partial applications currently?</b>	Yes	Yes	Yes
<b>Other unique features</b>	Real Color Scans, HD Photo, Scans Hard Palate, 3 second live bite, shade matching, Direct integration with Dextrix powered by DDX O/S.	Small tip, single button control, 280-gram weight, automatic and manual margin line creations, HD camera, scan replay.	Scans in Color & Can operate from Surface Pro Tablet; Can integrate with any computer in the office that has minimum specs. Direct integration to DPMS, means their scans can be reviewed or sent from ANY COMPUTER in the office.
<b>Compatible with implant scan bodies for custom abutments?</b>	Yes, 3i Encode and all others	Yes, utilizes AI to automatically match and align scan abutments and scan bodies with scan data.	Yes, 3i Encode and all others

Company	Align Technology	Dentsply Sirona
<b>Product</b>	<b>iTero® Element 5D</b> 	<b>CEREC® PRIMESCAN</b> 
<b>Scanner File Output</b>	STL	DXD, RST, 3D PDF, STL (via lab software)
<b>Scanner Type</b>	Parallel Confocal; Laser	Streaming Capture
<b>Scanner Accuracy (in microns)</b>	<9	Dentsply Sirona does not supply accuracy in microns
<b>Infection Control</b>	Disposable Lense Cap	High Level Disinfection & Dry Heat Sterilization
<b>Tooth Surface Spray Required?</b>	No	No
<b>Per Scan Click Fee?</b>	No	No
<b>Subscription (data package and/or software updates)</b>	Starting in yr 2, \$360/month (\$4,320 annually) all inclusive monthly subscription	No subscription fees
<b>Lapsed Time (min) until digital model file is available</b>	15 to 30 mins. Varies depending on office's internet connection	0 to 5; Varies Depending on office's internet connection
<b>Hardware Warranty?</b>	Included with subscription	1 year Included with purchase. Optional \$199/month extends warranty up to 7 yrs coverage for hardware warranty and design software updates.
<b>In-House Training</b>	No in-office training. 6-8 hours of virtual classroom style training. In-office training can be purchased.	Training included with intial purchase price. No obligation - additional training is offered
<b>Suggested Retail Price</b>	\$29,999	\$30,000
<b>Distributors</b>	Direct and Patterson Dental	Patterson Dental and Henry Schein
<b>Ortho applications currently?</b>	Invisalign, OrthoCAD; SureSmile	Yes-MTM Clear Aligners, Clear Correct, and Invisalign
<b>Partial applications currently?</b>	Yes	Yes
<b>Other unique features</b>	CDT inspects impression file and marks margins before manufacturing lab receives them.	Full-color scanning, small camera tip, restorative and ortho scanning solutions, direct integration into Dolphin Imaging. 3 machine configurations for maximum flexibility and portability.
<b>Compatible with implant scan bodies for custom abutments?</b>	Yes, 3i Encode and all others	Yes-most major implant companies utilizing Dentsply Sirona TiBase/ScanPost and Atlantis workflows (2017)